


Neutron diffraction and imaging on battery systems


senyshyn@frm2.tum.de

MLZ is a cooperation between:

Li-ion battery: principle of operation




Market price \$/kWh based on Li-ion technology





Courtesy: Gleb Yushin

Projected production of Li-ion batteries in Europe



Time evolution of 18650 cell capacity





Moore's Law does not hold true for battery development



- 120 mAh/annum average increase rate over almost 25 years
- ~ 3.5 %/year relative increase
- Since 2012 the capacity increase is achieved by voltage increase and introduction of Si to graphite anodes

Materials for battery applications


Anode materials



Cathode materials




Lithium electrolytes




x=0.50

z=1.00



Different mechanisms of Li-ion battery degradation



Copyright Christoph Birkl 2014

Neutron-based experimental techniques with proven relevance\impact in battery research

Neutron diffraction: detail of crystal structure, localisation and quantification of lithium; microstructural studies; phase analysis.

Neutron imaging: lithium distribution, gas formation, electrolyte dynamics;

Small-angle neutron scattering: in-situ materials morphology and fracturing upon cell fatigue;

Quasielastic neutron scattering: in-situ structure and mobility of electrolytes in Li-ion batteries;

Reflectometry: studies of solid-electrolyte interphase; studies of lithiation in amorphous silicon; solid-liquid interfaces;


Neutron depth profiling: nanometer sensitive probe of lithium concentration in electrode materials;

Positron spectroscopy: charge- and fatigue-induced defect formation;


Neutron and Prompt gamma activation analysis: non-destructive and simultaneous elemental/isotope analysis;

Gas evolution in pouch cells studied by neutron radiography

LFP\|LTO cell




LFP\|graphite cell



B. Michalak et al., *Sci. Rep.* 5 (2015) 15627

4D imaging on lithium-batteries

Pristine/0 s Partly dc/1500 s
-225.71 mAh -580.55 mAh





R.F. Ziesche et al., Nature Communications 11 (2020) 777

Simultaneous neutron radiography and diffraction data collection on 18650-type cell cycled up-side-down

Fresh cell


Fatigued cell



https://www.youtube.com/watch?v=ICPzHO_1nQ8




Why graphite?

Neutron, $\lambda=1.5482 \text{ \AA}$



High-energy synchrotron,



$\lambda=0.20708 \text{ \AA}$



Simultaneous neutron radiography and diffraction data collection on 18650-type cell cycled up-side-down




Fresh cell

Fatigued cell




https://www.youtube.com/watch?v=ICPzHO_1nQ8



Spatially-resolved neutron diffraction and current distribution in Li-ion batteries



Selected diffraction patterns






Lithium distribution in the middle of 18650-type cell



Spatially-resolved TOF neutron diffraction

OAK RIDGE
National Laboratory




Fresh

120 cycles

200 cycles

400 cycles



D. Petz et al., J. Power Sources 448 (2020) 227466

Spatially-resolved diffraction using conical slits


Gauge volume

Neutron diffraction





Gauge volume

Diffraction tomography





Principle of diffraction tomography





Lithium distribution in the graphite anode of 18650-type lithium ion battery

Spatially resolved neutron diffraction



Neutron diffraction tomography*



Summary

- Perspectives for neutron powder diffraction and spatially-resolved diffraction

More flux, better resolution



- Neutron diffraction tomography

Low-divergent \parallel monochromatic neutron beam with submillimeter focusing

- Neutron imaging in battery research

More flux, better resolution, higher neutron energy

pinhole neutron diffraction



Acknowledgment

Dr. M.J. Mühlbauer (DPMA), Dr. O. Dolotko (Ames Lab)

Dr. V. Baran (DESY), D. Petz (TUM), V. Kochetov (Uni Rostock)

STRESS-SPEC team: Dr. M. Hoffmann (FRM II, TUM), Dr. W. Gan (HZG)

Dr. J. Rebell-Kornmeier, K. Braun

ILL: Dr. T. Pirling (SALSA, ILL)

KIT: Dr. M. Heere, Dr. M. Knapp, Prof. H. Ehrenberg

Berkeley Lab: Dr. R. Kostecki



Bundesministerium
für Bildung
und Forschung


DFG Deutsche
Forschungsgemeinschaft




FRM II
Forschungs-Neutronenquelle
Heinz Maier-Leibnitz

Neutron CT reconstruction from 18650-type LCO|C cell


ANTARES@FRMII



3.0 V




Aluminum for Cathode




Anode Electrodes



4.2 V



Neutron CT reconstruction from 18650-type LCO|C cell





SOC

— SOC 0
— SOC 100
— difference




Fatigue

— SOC 100, fresh
— SOC 100, fatigued
— difference



Evolution of the neutron diffraction signal (background subtracted) upon cooling of LP30 electrolyte filled in a thin-wall vanadium container.

Chosen electrolyte: LP30, EC+DMC+1M LiPF₆; Melting temperature: ca.250 K



Distribution of lithium and electrolyte concentration in fresh and aged 18650-type cells

